

Shoprider VENICE

ASK FOR THE BRAND **SHOPRIDER**



SHOPRIDER[®]
www.shoprider.com.au
mobility

RELIABLE MOBILITY GUARANTEED

Shoprider

VENICE

SHOPRIDER®
www.shoprider.com.au *mobility*

The SHOPRIDER™ **VENICE** powerchair features World Leading Battery Technology to reduce the overall weight of the unit. The **Venice** is a versatile powerchair ideally suited to the needs of the portable user. The **Venice** can be dismantled for transport by using the quick disconnect lever in seconds. The **Venice** features a comfortable mid back seat, adjustable height armrests and when combined with the powerful motors it will deliver outstanding safety & comfort.

SPECIFICATIONS

Length	725mm
Width	545mm
Height	1180mm
Battery Capacity -Lithium Iron Phosphate	10 Amp
#Maximum Range	up to 16km
Maximum Load Single User*	113Kgs
Charger on board	2Amp
Motor Size	500watt
Maximum Speed	6Kph
*Safe Climbing Angle	10 Degrees
*Max Climbing Angle	10 Degrees
Turning Radius	400mm
Wheel Size:	
- Front Caster	2 x 155mm
- Rear Caster	2 x 200mm

FEATURES

Lithium Iron Phosphate Battery
Mid Back Captain Seat
Colour Change Body Covers
Variable seat height
Adjustable arm rests
Controller interchangeable
Free wheel device
Easy Dismantle for Transport

Driver weight may exceed weight of unit.
Reduce speed when turning.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

VENICE



*(Single user)

10° (17%) - 90kg
6° (10%) - 100kg
0° (0%) - 113kg

SHOPRIDER®
www.shoprider.com.au *mobility*

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au

Your Authorised dealer
Coral Coast Mobility
Ph: 0899 215 954
Mob: 0447 215 954